



Quantitative Labeling Using Tandem Mass Tags


Rosa Viner
Omics marketing group
September 13, 2018



Introduction: TMT based relative quantification


Moving Beyond Qualitative Proteomics

Problem: Quantitative information about expression level of a protein is essential to understanding its biological role in response to change or disease.




Add another dimension to any experiment by determining the relative abundance of each identified protein

Isobaric Labeling/Tandem Mass Tags™ (TMT)*




* Ting, L. et al. 2011. Nature Methods 8: 937-940
Tandem Mass Tag and TMT are trademarks of Proteome Sciences plc.

How Does Isobaric Mass Tagging (TMT&iTRAQ) Work?




The tagged peptides behave exactly the same, except during fragmentation.





This spectra indicates that this protein is upregulated in patient 1 approximately 2 fold

Tandem Mass Tag and TMT are trademarks of Proteome Sciences plc.

TMT Labels are Indistinguishable




The Difference Only Appears in MSn




Synchronous Precursor Selection (SPS) for Accurate Quantification


Precursor Ion



CID MS², Ion Trap




HCD MS³, Orbitrap



- Available on the Thermo Scientific™ Orbitrap Fusion™ MS and Thermo Scientific™ Orbitrap Fusion™ Lumos™ MS
- Select multiple MS² precursors using a single fill and notched isolation waveform
- Improves the ratio accuracy and at the same time dramatically boosts sensitivity

Multinotch MS³ Quantification is Accurate and Sensitive






Synchronous Precursor Selection

Number of Precursors: 10

MS Isolation Window (m/z): 2


MS2 Isolation Window (m/z): 2



* TMT MS³ Peptides-Quan template in Method Editor

TMT11plex


Reporter Ion
Mass







- TMT11-131C can be used in combination with TMT10plex reagents to multiplex 11 different samples for MS analysis
- 11plex data analysis is supported by Proteome Discoverer 2.1-2.3



High Resolving Power is Required for Accurate Quantification of the TMT11 plex



TMT Multiplexing Workflow for Precise Data in Less Time



Sample Labeling and Preparation

LC-MS/MS (SPS MS3) Analysis


Data Analysis



Thermo Scientific™ TMT™
11-Plex Reagents



Thermo Scientific™ Orbitrap Fusion™
Lumos™ Tribrid™ MS with Method Templates



Thermo Scientific™ Proteome
Discoverer™ Software and
ProteinCenter™ Software

Unique workflow with potential for massive throughput – *improvements in instrument technology*



ThermoFisher
S C I E N T I F I C

TMT QC assay - TMT 11 yeast triple knock out standard

The world leader in serving science

- MS and LC method optimization
- QC of mass spec and LC
- PD analysis optimization

“TMT standard will not solve the interference problem, it can accurately and sensitively measure its effects”

J.Paulo et al, JASMAS, 2016, 1620-1625



© American Society for Mass Spectrometry, 2016




J. Am. Soc. Mass Spectrom. (2016)
DOI: 10.1007/s13361-016-1434-9

RESEARCH ARTICLE

A Triple Knockout (TKO) Proteomics Standard for Diagnosing Ion Interference in Isobaric Labeling Experiments

Joao A. Paulo, Jeremy D. O'Connell, Steven P. Gygi

Department of Cell Biology, Harvard Medical School, Boston, MA 02115, USA







Abstract. Isobaric labeling is a powerful strategy for quantitative mass spectrometry-based proteomic investigations. A complication of such analyses has been the co-isolation of multiple analytes of similar mass-to-charge resulting in the distortion of relative protein abundance measurements across samples. When properly implemented, synchronous precursor selection and triple-stage mass spectrometry (SPS-MS3) can reduce the occurrence of this phenomenon, referred to as ion interference. However, no diagnostic tool is available currently to rapidly and accurately assess ion interference. To address this need, we developed a multiplexed tandem mass tag (TMT)-based standard, termed the triple knockout (TKO). This standard is comprised of three yeast proteomes in triplicate, each from a strain deficient in a highly abundant protein (Met6, Pfk2, or Ura2). The relative abundance patterns of these proteins, which can be inferred from dozens of peptide measurements can demonstrate ion interference in peptide quantification. We expect no signal in channels where the protein is knocked out, permitting maximum sensitivity for measurements of ion interference against a null background. Here, we emphasize the need to investigate further ion interference-generated ratio distortion and promote the TKO standard as a tool to investigate such issues.

Keywords: MS standard, MultiNotch, TMT, Orbitrap Fusion, Lumos, Ion interference, SPS-MS3

Received: 12 April 2016/Revised: 30 May 2016/Accepted: 31 May 2016





Difference between Proteomics TKO Standard and Pierce TMT TKO standard

- TKO standard (Paulo et al, 2016)
 - Paulo's TKO standard uses met6(protein rank19), pfk2(114) and ura2(244) yeast strains.
 - Paulo et al uses TMT9 to label peptides from those three strains.



- Pierce TMT11 TKO Standard

- Pierce TMT11 standard uses met6,his4(213),ura2(244), and BY4742 parental yeast strains.
- The standard uses BY4742 parental strain labeled with 131N and 131C as control channels.





ThermoFisher
S C I E N T I F I C

Experimental Set up

The world leader in serving science

TMT11 TKO Sample Reconstitution

- Take vial out of the freezer and bring to the room temperature, 15 min
- Add 40 ul of 0.1%TFA/5% AcN in Water, 500 ng/ul
- Or add 40 ul of 5% DMSO/1% FA in Water, 500 ng/ul
- Incubate at RT for 15 min and transfer to autosampler vial
- Don't keep more than 1 week at 4 C

LC Method: EASY-Spray™ C18 50cm column & UltiMate™ 3000 RSLCnano UHPLC

50 min gradient

| No | Time | Flow [µl/min] | %B | Curve |
|----|----------------|---------------|------|-------|
| 1 | 0.000 | Equilibration | | |
| 2 | 0.000 | 0.300 | 2.0 | 5 |
| 3 | <i>New Row</i> | | | |
| 4 | 0.000 | Run | | |
| 5 | 14.000 | 0.300 | 2.0 | 5 |
| 6 | 17.000 | 0.300 | 4.0 | 5 |
| 7 | 67.000 | 0.300 | 28.0 | 5 |
| 8 | 70.000 | 0.300 | 65.0 | 5 |
| 9 | 75.000 | 0.300 | 65.0 | 5 |
| 10 | 77.000 | 0.300 | 4.0 | 5 |
| 11 | <i>New Row</i> | | | |
| 12 | 100.000 | Stop Run | | |

120 min gradient

| No | Time | Flow [µl/min] | %B | Curve |
|----|----------------|---------------|------|-------|
| 1 | 0.000 | Equilibration | | |
| 2 | 0.000 | 0.300 | 2.0 | 5 |
| 3 | <i>New Row</i> | | | |
| 4 | 0.000 | Run | | |
| 5 | 14.000 | 0.300 | 2.0 | 5 |
| 6 | 17.000 | 0.300 | 4.0 | 5 |
| 7 | 100.000 | 0.300 | 16.0 | 5 |
| 8 | 145.000 | 0.300 | 25.0 | 5 |
| 9 | 150.000 | 0.300 | 65.0 | 5 |
| 10 | 158.000 | 0.300 | 65.0 | 5 |
| 11 | 160.000 | 0.300 | 4.0 | 5 |
| 12 | <i>New Row</i> | | | |
| 13 | 185.000 | Stop Run | | |



Solvent A: 0.1% formic acid
 Solvent B: 100 % Acetonitrile, 0,1% formic acid
 Flow rate: 300 nL/min
 Injection volume: 1-2 µL

Direct injection



LC Method: EASY-Spray™ C18 50cm column & UltiMate™ 3000 RSLCnano UHPLC

Loading pump conditions

| Time(min) | Flow, ul/min |
|-----------|--------------|
| 0 | 20 |
| 100, 185 | 20 |

50 min gradient

| No | Time | Flow [µl/min] | %B | Curve |
|----|----------------|---------------|------|-------|
| 1 | 0.000 | Equilibration | | |
| 2 | 0.000 | 0.300 | 2.0 | 5 |
| 3 | <i>New Row</i> | | | |
| 4 | 0.000 | Run | | |
| 5 | 14.000 | 0.300 | 2.0 | 5 |
| 6 | 17.000 | 0.300 | 4.0 | 5 |
| 7 | 67.000 | 0.300 | 28.0 | 5 |
| 8 | 70.000 | 0.300 | 65.0 | 5 |
| 9 | 75.000 | 0.300 | 85.0 | 5 |
| 10 | 77.000 | 0.300 | 4.0 | 5 |
| 11 | <i>New Row</i> | | | |
| 12 | 100.000 | Stop Run | | |

120 min gradient

| No | Time | Flow [µl/min] | %B | Curve |
|----|----------------|---------------|------|-------|
| 1 | 0.000 | Equilibration | | |
| 2 | 0.000 | 0.300 | 2.0 | 5 |
| 3 | <i>New Row</i> | | | |
| 4 | 0.000 | Run | | |
| 5 | 14.000 | 0.300 | 2.0 | 5 |
| 6 | 17.000 | 0.300 | 4.0 | 5 |
| 7 | 100.000 | 0.300 | 16.0 | 5 |
| 8 | 145.000 | 0.300 | 25.0 | 5 |
| 9 | 150.000 | 0.300 | 85.0 | 5 |
| 10 | 158.000 | 0.300 | 85.0 | 5 |
| 11 | 160.000 | 0.300 | 4.0 | 5 |
| 12 | <i>New Row</i> | | | |
| 13 | 185.000 | Stop Run | | |



Solvent A: 0.1% formic acid
 Solvent B: 100 % Acetonitrile, 0,1% formic acid
 Flow rate: 300 nL/min
 Injection volume: 1-2 µL

Trap loading



LC Method using EASY-Spray™ C18 50cm column & EASY-nLC™ 1200 HPLC

50 min

120 min

| Time (min) | Flow(nL/min) | %B | Time (min) | Flow(nL/min) | %B |
|------------|--------------|----|------------|--------------|----|
| 0 | 300 | 5 | 0 | 300 | 5 |
| 5 | 300 | 10 | 5 | 300 | 8 |
| 55 | 300 | 40 | 125 | 300 | 40 |
| 60 | 300 | 90 | 130 | 300 | 90 |
| 70 | 300 | 90 | 140 | 300 | 90 |



Solvent A: 0.1% formic acid
 Solvent B: 80 % Acetonitrile, 0,1% formic acid
 Flow rate: 300 nL/min
 Injection volume: 1-2 µL


LC Method using EASY-Spray™ C18 50cm column & EASY-nLC™ 1000 HPLC

50 min

120 min

| Time (min) | Flow(nL/min) | %B | Time (min) | Flow(nL/min) | %B |
|------------|--------------|----|------------|--------------|----|
| 0 | 300 | 5 | 0 | 300 | 5 |
| 5 | 300 | 7 | 5 | 300 | 7 |
| 55 | 300 | 32 | 125 | 300 | 32 |
| 60 | 300 | 90 | 130 | 300 | 90 |
| 70 | 300 | 90 | 140 | 300 | 90 |

Solvent A: 0.1% formic acid
 Solvent B: 100 % Acetonitrile, 0,1% formic acid
 Flow rate: 300 nL/min
 Injection volume: 1 µL




Instrument Method Settings: QExactive Classic and QExactive Plus

| Properties | QE classic 50 min | QE classic 120 min | QE + 50 min | QE + 120 min |
|------------------------|-------------------|--------------------|-------------|--------------|
| Resolution Full MS | 70000 | 70000 | 70000 | 70000 |
| AGC target Full MS | 3e6 | 3e6 | 3e6 | 3e6 |
| MS max IT, ms | 50 | 50 | 50 | 50 |
| Scan range, <i>m/z</i> | 350-1500 | 350-1500 | 350-1500 | 350-1500 |
| Loop count | 15 | 15 | 15 | 15 |
| MS2 resolution | 35000 | 35000 | 35000 | 35000 |
| MS2 AGC target | 1e5 | 1e5 | 1e5 | 1e5 |
| MS2 max IT, ms | 120 ms | 250 ms | 100 ms | 120 ms |
| Isolation Window , Th | 1.2 | 1.2 | 0.7 | 0.7 |
| NCE, % | 32-34 | 32-34 | 32-34 | 32-34 |
| Intensity threshold | 1e4 | 1e4 | 1e4 | 1e4 |
| Peptide match | preferred | preferred | preferred | preferred |
| Dynamic exclusion, s | 20 s | 45 s | 20 s | 30 s |
| First mass, <i>m/z</i> | 110 | 110 | 110 | 110 |




Instrument Method Settings: QExactive HF and QExactive HF X

| Properties | QE HF 50 min | QE HF 120 min | QE HF X 50 min | QE HF X 120 min |
|----------------------|-------------------|-------------------|---------------------------|---------------------------|
| Resolution Full MS | 120000 | 120,000 | 120000 | 120000 |
| AGC target Full MS | 3e6 | 3e6 | 3e6 | 3e6 |
| MS max IT, ms | 50 | 50 | 50 | 50 |
| Scan range, m/z | 350-1500 | 350-1500 | 350-1500 | 350-1500 |
| Loop count | 20 | 15 | 20 | 15 |
| MS2 resolution | 60000(2.9-45000) | 60000(2.9- 45000) | 45000 | 45000 |
| MS2 target | 1e5 | 1e5 | 1e5 | 1e5 |
| MS2 max IT, ms | 96 | 120 | 86 | 96 |
| Isolation Window, Th | 0.4 m/z | 0.7 m/z | 0.7 m/z | 0.7 m/z |
| NCE, % | 32-34 | 32-34 | 32-34 | 32-34 |
| Intensity threshold | 1e4 | 1e4 | 1e4 | 1e4 |
| Peptide match | preferred | preferred | preferred , single charge | preferred , single charge |
| Dynamic exclusion, s | 20 | 30 | 20 | 30 |
| First mass, m/z | 110 | 110 | 110 | 110 |




Instrument Method Settings: Fusion, Tune3.1 Templates are Available in Method Editor

| Properties | Fusion SPS 50 min | Fusion SPS 120 min | Fusion MS2 50 min | Fusion MS2 120 min |
|--|-------------------------------|-------------------------------|------------------------|--------------------|
| Resolution Full MS | 120000 | 120000 | 120000 | 120000 |
| AGC target Full MS | 4e5 | 4e5 | 4e5 | 4e5 |
| MS max IT, ms | 50 | 50 | 50 | 50 |
| Scan range, m/z | 375-1500 | 375-1500 | 375-1500 | 375-1500 |
| Top Speed, s | 2 | 3 | 2 | 3 |
| MS2 max IT, ms | 50 | 50 | 105 | 120 |
| MS2 Isolation window, Th | 1.2(2)-0.7(3)-0.5 (4+) | 1.2(2)-0.7(3)-0.5 (4+) | 1.2(2)-0.7(3)-0.5 (4+) | 0.7(2-3)-0.5 (4+) |
| MS2 NCE, % | 35 | 35 | 38-40 | 38-40 |
| MS2 Intensity threshold | 5e3 | 5e3 | 5e4 | 5e4 |
| Dynamic exclusion, s | 45, single charge | 60, single charge | 45, single charge | 60, single charge |
| MS2 Resolution | turbo | turbo | 50000 | 50000 |
| MS2 AGC target | 1e4 | 1e4 | 1e5 | 1e5 |
| MS3 AGC target | 1e5 | 1e5 | | |
| SPS Isolation window, Th | 1.3(2)-0.7(3)-0.5 (4+) | 1.3(2)-0.7(3)-0.5 (4+) | | |
| SPS NCE, % | 65 | 65 | | |
| SPS max IT, ms | 105 | 120 | | |
| SPS settings: # notches, mass range | 5-10-10 <i>m/z</i> 110-500 | 5-10-10 <i>m/z</i> 110-500 | <i>m/z</i> 110 | <i>m/z</i> 110 |



Instrument Method Settings: Lumos, Tune3.1 Templates are Available in Method Editor

| Properties | Fusion SPS 50 min | Fusion SPS 120 min | Fusion MS2 50 min | Fusion MS2 120 min |
|--|-------------------------------|-------------------------------|------------------------|--------------------|
| Resolution Full MS | 120000 | 120000 | 120000 | 120000 |
| AGC target Full MS | 4e5 | 4e5 | 4e5 | 4e5 |
| MS max IT, ms | 50 | 50 | 50 | 50 |
| Scan range, m/z | 375-1500 | 375-1500 | 375-1500 | 375-1500 |
| Top Speed, s | 2 | 3 | 2 | 3 |
| MS2 max IT, ms | 50 | 50 | 86 | 96 |
| MS2 Isolation window, Th | 1.2(2)-0.7(3)-0.5 (4+) | 1.2(2)-0.7(3)-0.5 (4+) | 1.2(2)-0.7(3)-0.5 (4+) | 0.7(2-3)-0.5 (4+) |
| MS2 NCE, % | 35 | 35 | 38-40 | 38-40 |
| MS2 Intensity threshold | 5e3 | 5e3 | 5e4 | 5e4 |
| Dynamic exclusion, s | 45, single charge | 60, single charge | 45, single charge | 60, single charge |
| MS2 Resolution | turbo | turbo | 50000 | 50000 |
| MS2 AGC target | 1e4 | 1e4 | 1e5 | 1e5 |
| MS3 AGC target | 1e5 | 1e5 | | |
| SPS Isolation window, Th | 1.3(2)-0.7(3)-0.5 (4+) | 1.3(2)-0.7(3)-0.5 (4+) | | |
| SPS NCE, % | 65 | 65 | | |
| SPS max IT, ms | 86 | 105 | | |
| SPS settings: # notches, mass range | 5-10-10 <i>m/z</i> 110-500 | 5-10-10 <i>m/z</i> 110-500 | <i>m/z</i> 110 | <i>m/z</i> 110 |







ThermoFisher
S C I E N T I F I C

PD 2.2 workflow set up

I. Quantification Method Add Lot Specific Correction Factors




1



2

| Mass Tag | Reporter Ion | -2 | -1 | Monoisotopic | +1 | +2 |
|-------------------------|--------------|------|------|--------------|------|------|
| TMT ¹⁰ -126 | 126.127726 | 0% | 0% | 100% | 8.5% | 0.3% |
| TMT ¹⁰ -127N | 127.124761 | 0% | 0.5% | 100% | 8.5% | 0.3% |
| TMT ¹⁰ -127C | 127.131081 | 0% | 0.5% | 100% | 7.3% | 0.4% |
| TMT ¹⁰ -128N | 128.128116 | 0% | 0.6% | 100% | 7.2% | 0.2% |
| TMT ¹⁰ -128C | 128.134436 | 0% | 1.3% | 100% | 6.3% | 0.2% |
| TMT ¹⁰ -129N | 129.131471 | 0% | 1.6% | 100% | 6.2% | 0% |
| TMT ¹⁰ -129C | 129.137790 | 0% | 2.5% | 100% | 5.0% | 0.1% |
| TMT ¹⁰ -130N | 130.134825 | 0% | 2.7% | 100% | 5.0% | 0% |
| TMT ¹⁰ -130C | 130.141145 | 0.1% | 2.9% | 100% | 4.0% | 0% |
| TMT ¹⁰ -131 | 131.138180 | 0.1% | 3.1% | 100% | 3.9% | 0% |

3



4

| Mass Tag | Reporter Ion Mass | - 2 | - 1 | Main | + 1 | + 2 | Active |
|----------|-------------------|-----|-----|------|-----|-----|-------------------------------------|
| 126 | 126.127726 | 0 | 0 | 100 | 8.5 | 0.3 | <input checked="" type="checkbox"/> |
| 127N | 127.124761 | 0 | 0.5 | 100 | 8.5 | 0.3 | <input checked="" type="checkbox"/> |
| 127C | 127.131081 | 0 | 0.5 | 100 | 7.3 | 0.4 | <input checked="" type="checkbox"/> |
| 128N | 128.128116 | 0 | 0.6 | 100 | 7.2 | 0.2 | <input checked="" type="checkbox"/> |
| 128C | 128.134436 | 0 | 1.3 | 100 | 6.3 | 0.2 | <input checked="" type="checkbox"/> |
| 129N | 129.131471 | 0 | 1.6 | 100 | 6.2 | 0 | <input checked="" type="checkbox"/> |
| 129C | 129.13779 | 0 | 2.5 | 100 | 5 | 0.1 | <input checked="" type="checkbox"/> |
| 130N | 130.134825 | 0 | 2.7 | 100 | 5 | 0 | <input checked="" type="checkbox"/> |
| 130C | 130.141145 | 0 | 2.9 | 100 | 4 | 0 | <input checked="" type="checkbox"/> |
| 131N | 131.13818 | 0 | 3.1 | 100 | 3.9 | 0 | <input checked="" type="checkbox"/> |
| 131C | 131.144499 | 0 | 1.4 | 100 | 2.9 | 0 | <input checked="" type="checkbox"/> |

TMT: Main peaks are always 100%

OK Cancel Help

1. Select "Maintain Quantification Methods"
2. Create new method using TMT 11plex template
3. Add Correction factors from Product data sheet
4. Save new method as TMT11TKOlotXXX standard

II. Study Set up

Thermo Proteome Discoverer 2.2.0.388

File View Administration Tools Window Help

Start Page x Study: test x Administration x F1_20170817_FL_HeLa_lug_OT_120min_high_charge_BIH x F1_20170817_FL_HeLa_lug_OT_120min_low_charge_BIH x hela500ng6lg2 x

Proteome Discoverer 2.2

Start

New Study/Analysis... 1

Open Study...

Open Result...

Recent Studies Recent Results 7

New Study and Analysis

Study Name: TMT 11 TKO standard

Study Root Directory: E:\TMOnevstandard 3

Processing Workflow: uan_SPS_MS3_SequestHT_Percolator.pdProcessingWF

4

Consensus Workflow: (empty workflow)

5

- (empty workflow)
- CWF_BasicXlink.pdConsensusWF
- PMI-Byonic Template.pdConsensusWF
- ConsensusWF \ CWF_Basic.pdConsensusWF
- ConsensusWF \ CWF_Basic_Annotation.pdConsensusWF
- ConsensusWF \ CWF_Comprehensive_Enhanced Annotation.pdConsensusWF
- ConsensusWF \ CWF_Comprehensive_Enhanced Annotation_LFQ_and_Precursor_Quan.pdConsensusWF
- ConsensusWF \ CWF_Comprehensive_Enhanced Annotation_Quan.pdConsensusWF
- ConsensusWF \ CWF_Comprehensive_Enhanced Annotation_Quan_Results export.pdConsensusWF
- ConsensusWF \ CWF_Comprehensive_Enhanced Annotation_Reporter_Quan.pdConsensusWF
- ProSightPD 1.1 for PD 2.2 and PSpC 4.0 Templates \ ProSightPD Bottom Up.pdConsensusWF
- ProSightPD 1.1 for PD 2.2 and PSpC 4.0 Templates \ ProSightPD HI HI.pdConsensusWF
- ProSightPD 1.1 for PD 2.2 and PSpC 4.0 Templates \ ProSightPD LO HI.pdConsensusWF
- ProSightPD 1.1 for PD 2.2 and PSpC 4.0 Templates \ ProSightPD MED HI.pdConsensusWF

1. Select New Study
2. Create Study Name
3. Select Study Directory
4. Select Processing workflow
5. Select Consensus workflow
6. Select Quan.method and control channel
7. Add files

Quantification Method: TMT11Universal

Select Control Channel:

| | |
|-------------------------------|--|
| <input type="checkbox"/> 126 | <input type="checkbox"/> 129C |
| <input type="checkbox"/> 127N | <input type="checkbox"/> 130N |
| <input type="checkbox"/> 127C | <input type="checkbox"/> 130C |
| <input type="checkbox"/> 128N | <input type="checkbox"/> 131N |
| <input type="checkbox"/> 128C | <input checked="" type="checkbox"/> 131C |
| <input type="checkbox"/> 129N | |

6

II. Study Set up

Analysis Configuration Window

Quantification Methods:

- Resolution TMTe 6plex
- 2plex (Arg10, Lys6)
- 2plex (Arg10, Lys8)
- 2plex (Ile6)
- SILAC 3plex (Arg6, Lys4 | Arg10, Lys8)
- SILAC 3plex (Arg6, Lys6 | Arg10, Lys8)
- TMT 10plex
- TMT 11plex
- TMT 2plex
- TMT11Universal**

Study Factors:

- Biological Replicate Factor
- Categorical Factor** (1)
- Numerical Factor

Yeast Strain Configuration

Yeast Strain: his4 met6 parental ura2 (2)

1. Select categorical factor
2. Create study factors
3. Set study factors, controls for each of the quan.channels per file

Thermo Proteome Discoverer 2.2.0.388

| ID | Name | File Type | Quan Method | Sample Information |
|----|-------------------|-----------|----------------|--|
| F1 | TKOTT11_1ms2_1 | .raw | TMT11Universal | Sample Type: [Sample], Yeast Strain: [n/a] |
| F2 | TKOTT11_1ms2_1_42 | .raw | TMT11Universal | Sample Type: [Sample], Yeast Strain: [n/a] |
| F3 | TKOTT11_1ms2_2 | .raw | TMT11Universal | Sample Type: [Sample], Yeast Strain: [n/a] |
| F4 | TKOTT11_1ms3_1 | .raw | TMT11Universal | Sample Type: [Sample, Control], Yeast Strain: [met6, his4, ura2, parental] |

| Sample | Sample Identifier | Sample Type | Quan Channel | Yeast Strain |
|--------|-------------------------|-------------|--------------|--------------|
| 4 | TKOTT11_1ms3_1 - [126] | Sample | 126 | met6 |
| 46 | TKOTT11_1ms3_1 - [127N] | Sample | 127N | met6 |
| 47 | TKOTT11_1ms3_1 - [127C] | Sample | 127C | met6 |
| 48 | TKOTT11_1ms3_1 - [128N] | Sample | 128N | his4 |
| 49 | TKOTT11_1ms3_1 - [128C] | Sample | 128C | his4 |
| 50 | TKOTT11_1ms3_1 - [129N] | Sample | 129N | his4 |
| 51 | TKOTT11_1ms3_1 - [129C] | Sample | 129C | ura2 |
| 52 | TKOTT11_1ms3_1 - [130N] | Sample | 130N | ura2 |
| 53 | TKOTT11_1ms3_1 - [130C] | Sample | 130C | ura2 |
| 54 | TKOTT11_1ms3_1 - [131N] | Control | 131N | parental |
| 55 | TKOTT11_1ms3_1 - [131C] | Control | 131C | parental |

III. Search Parameters For Processing Workflow

Download SwissProt yeast database, taxonomy ID 4932

Parameters of 'Spectrum Files RC'

Show Advanced Parameters

1. Search Settings

Protein Database Saccharomyces cerevisiae S288C (SwissProt)

Enzyme Name Trypsin (Full)

1. Dynamic Modification None

Static Peptide N-Terminus TMT6plex / +229.163 Da (Any N-Terminus)

1. Static Modification TMT6plex / +229.163 Da (K)

Weight of c ions 0

Weight of x ions 0

Weight of y ions 1

Weight of z ions 0

4. Dynamic Modifications

Max. Equal Modifications 3

1. Dynamic Modification Oxidation / +15.995 Da (M)

2. Dynamic Modification None

3. Dynamic Modification None

4. Dynamic Modification None

5. Dynamic Modification None

6. Dynamic Modification None

5. Dynamic Modifications (peptide terminus)

1. N-Terminal Modification None

2. N-Terminal Modification None

3. N-Terminal Modification None

1. C-Terminal Modification None

2. C-Terminal Modification None

3. C-Terminal Modification None

6. Dynamic Modifications (protein terminus)

1. N-Terminal Modification Acetyl / +42.011 Da (N-Terminus)

2. N-Terminal Modification None

3. N-Terminal Modification None

1. C-Terminal Modification None

2. C-Terminal Modification None

3. C-Terminal Modification None

7. Static Modifications

Peptide N-Terminus TMT6plex / +229.163 Da (Any N-Terminus)

Peptide C-Terminus None

1. Static Modification Carbamidomethyl / +57.021 Da (C)


2. Static Modification TMT6plex / +229.163 Da (K)

3. Static Modification None

4. Static Modification None


5. Static Modification None

6. Static Modification None




Use 0.02 Da tolerance for MS2 methods and 1.2 Da for SPS methods

III. Consensus Workflow



- 1. **General Quantification Settings**
 - Peptides to Use Unique
 - Consider Protein Group True
 - Reject Quant Results False
- 2. **Reporter Quantification**
 - Reporter Abundance Automatic
 - Apply Quant Value Cor False
 - Co-Isolation Threshold 50 **75 for MS3**
 - Average Reporter S/N 10
- 3. **Normalization and Scaling**
 - Normalization Mode None
 - Proteins For Normalization
 - Scaling Mode On Controls Average
- 4. **Exclude Peptides from Protein Quantification**
 - 1. Excluded Peptide None
 - 2. Excluded Peptide None
 - 3. Excluded Peptide None
 - N-Terminal Excluded None
- 5. **Quan Rollup_Hypothesis Testing**
 - Ratio Calculation Summed Abundance Based
 - Maximum Allowed Fold 100
 - Imputation Mode None
 - Hypothesis Test ANOVA (Individual Proteins)
- 6. **Quan Ratio Distributions**
 - 1st Fold Change Thre 2
 - 2nd Fold Change Thre 4
 - 3rd Fold Change Thre 6

III. Consensus Workflow: New In PD 2.3 SPS Mass Matches



- 1. **General Quantification Settings**
 - Peptides to Use Unique
 - Consider Protein Gr: True
 - Reject Quan Result: False
- 2. **Reporter Quantification**
 - Reporter Abundance: Automatic
 - Apply Quan Value C: True
 - Co-Isolation Thresh: 75
 - Average Reporter S: 10
 - SPS Mass Matches: 65
- 3. **Normalization and Scaling**
 - Normalization Mode: None
 - Proteins For Normali:
 - Scaling Mode: On Controls Average
- 4. **Exclude Peptides from Protein Quantification**
 - For Normalization: Use All Peptides
 - For Protein Roll-Up: Use All Peptides
 - For Pairwise Ratios: Exclude Modified
 - 1. Considered Peptic: None
 - 2. Considered Peptic: None
 - 3. Considered Peptic: None
 - N-Terminal Consider: None
- 5. **Quan Rollup and Hypothesis Testing**
 - Protein Ratio Calcul: Protein Abundance Based
 - Maximum Allowed Fr: 100
 - Imputation Mode: None
 - Hypothesis Test: ANOVA (Individual Proteins)
- 6. **Quan Ratio Distributions**
 - 1st Fold Change Thr: 2
 - 2nd Fold Change Th: 4

SPS Mass Matches


Sequence: DGQIITPK, K9-TMT6plex (229.16293 Da), D1-TMT6plex (229.16293 Da)
 Charge: +2, Monoisotopic m/z: 715.93493 Da (+0.34 mmu/+0.48 ppm), MH+: 1430.86259 Da, RT: 72.4247 min,
 Identified with: Sequest HT (v1.17); XCorr:3.08, Percolator q-Value:1.2e-4, Percolator PEP:6.6e-4,
 Fragment match tolerance used for search: 0.6 Da

Fragment Matches

Value Type: Theo. Mass [Da]

Ion Series: Neutral Losses | Precursor Ions | Internal Fragments

| #1 | b ⁺ | b ²⁺ | Seq. | y ⁺ | y ²⁺ | #2 |
|----|----------------|-----------------|------------|----------------|-----------------|----|
| 1 | 345.19715 | 173.10221 | D-TMT6plex | | | 9 |
| 2 | 402.21862 | 201.61295 | G | 1086.67203 | 543.83965 | 8 |
| 3 | 530.27719 | 265.64223 | Q | 1029.65056 | 515.32892 | 7 |
| 4 | 643.36126 | 322.18427 | I | 901.59198 | 451.29963 | 6 |
| 5 | 744.40894 | 372.70811 | T | 788.50792 | 394.75760 | 5 |
| 6 | 857.49300 | 429.25014 | I | 687.46024 | 344.23376 | 4 |
| 7 | 958.54068 | 479.77398 | T | 574.37618 | 287.69173 | 3 |
| 8 | 1055.59344 | 528.30036 | P | 473.32850 | 237.16789 | 2 |
| 9 | | | K-TMT6plex | 376.27574 | 188.64151 | 1 |



III. Ratios Set up per Individual File

The screenshot displays the Thermo Proteome Discoverer 2.2.0.388 software interface. The main window is titled 'Study: TMT 11 TKO standard' and is currently in the 'Grouping & Quantification' tab. The interface is divided into several sections:

- Sample Group and Quan Ratio Specification:** This section contains 'Study Variables' with checkboxes for 'File', 'Quan Channel', 'Yeast Strain' (checked), and 'Sample Type'. Below this are 'Manual Ratio Generation' and 'Bulk Ratio Generation' sections, each with dropdown menus for 'Numerator' and 'Denominator' and an 'Add Ratio' button.
- Generated Sample Groups:** This section lists sample groups and their corresponding files. For example, the 'met6' group (F4) includes files 126, 127N, and 127C. The 'his4' group (F4) includes files 128N, 128C, and 129N. The 'ura2' group (F4) includes files 129C, 130N, and 130C.
- Generated Ratios:** This section shows a list of ratios, such as 'met6 / parental F4', 'his4 / parental F4', and 'ura2 / parental F4', each with a red 'X' icon to its left.
- Analysis Dialog Box:** A dialog box titled 'Analysis' is open on the right. It features a 'Consensus Step' section with a workflow 'CWF_Comprehensive_Enhanced Annotation_Quan' and a 'Processing Step' section with a workflow 'PW/F_QE_Reporter_Based_Quan_SequestHT_Percolator'. The 'Files for Analysis' section lists two files: 'F4 TKOTT11_1ms3_1 TMT11Universal Sample Type: [Sample, Control], Yeast S' and 'F5 TKOTT11_1ms3_2 TMT11Universal Sample Type: [Sample, Control], Yeast S'. The 'As Batch' checkbox is checked and circled in red.

III. Ratios Set up for Multiple files

Thermo Proteome Discoverer 2.2.0.388

File View Administration Tools Window Help

Start Page x Study: test x Administration x F1_20170817_FL_HeLa_lug_OT_120min_low_charge_BIH x Study: TMT 11 TKO standard x

Add Files Add Fractions Remove Files Open Containing Folder New Analysis Open Analysis Template

Study Definition Input Files Samples Analysis Results Workflows **Grouping & Quantification**

Sample Group and Quan Ratio Specification

Study Variables

- File
- Quan Channel
- Yeast Strain
- Sample Type

Manual Ratio Generation

Numerator:

Denominator:

Add Ratio

Bulk Ratio Generation

Denominators to be used:

- Yeast Strain : met6
- Yeast Strain : his4
- Yeast Strain : ura2
- Yeast Strain : parental

Add Ratios

Generated Sample Groups

met6

- 126 Sample met6 F4: TKOTT11_1ms3_1
- 127N Sample met6 F4: TKOTT11_1ms3_1
- 127C Sample met6 F4: TKOTT11_1ms3_1
- 126 Sample met6 F5: TKOTT11_1ms3_2
- 127N Sample met6 F5: TKOTT11_1ms3_2
- 127C Sample met6 F5: TKOTT11_1ms3_2

his4

- 128N Sample his4 F4: TKOTT11_1ms3_1
- 128C Sample his4 F4: TKOTT11_1ms3_1
- 129N Sample his4 F4: TKOTT11_1ms3_1
- 128N Sample his4 F5: TKOTT11_1ms3_2
- 128C Sample his4 F5: TKOTT11_1ms3_2
- 129N Sample his4 F5: TKOTT11_1ms3_2

Generated Ratios Clear All

- X met6 / parental
- X his4 / parental
- X ura2 / parental

Analysis As Batch Run Save X

Consensus Step X

Workflow: CWF_Comprehensive_Enhanced Annotation_Quan
Result File: TKOTT11_1ms3_1.pdResult

Child Steps: (1) Add

Processing Step Clone

Workflow: PWF_QE_Reporter_Based_Quan_SequestHT_Percolator
Result File: TKOTT11_1ms3_1.msf

Files for Analysis: (2) Clear All

- X F4 TKOTT11_1ms3_1 TMT11Universal Sample Type: [Sample, Control], Yeast S
- X F5 TKOTT11_1ms3_2 TMT11Universal Sample Type: [Sample, Control], Yeast S



ThermoFisher
S C I E N T I F I C

Results and Quality Control of LC-MS System

The world leader in serving science

Search and Quan Results: 1012 Protein Groups, 5777 Unique Peptides, 6372 PSMs

Thermo Proteome Discoverer 2.2.0.388

File View Administration Tools Window Help

Start Page x Study: TNO standard x TKOTT11_1ms3_2unique75 x TKOTT11_1ms3_2unique75-(1) x



Specialized Traces Result Statistics

| Protein Accession | Description | Coverage [%] | Sequence Coverage | # Peptides | # PSMs | Abundance Ratios | | | Abundances (Grouped) | | | | Abundances (Grouped) CV: [+] | | | | Found in Samples | | | |
|-------------------|--|--------------|-------------------|------------|--------|---------------------|---------------------|---------------------|----------------------|-------|-------|----------|--------------------------------|------|------|------|------------------|---------|---------|---------|
| | | | | | | (MetB) / (parental) | (His4) / (parental) | (Ura2) / (parental) | MetB | His4 | Ura2 | parental | 2.41 | 2.90 | 1.75 | 2.05 | 2.41 | 2.43 | 0.94 | 1.52 |
| P32861 | UTP-glucose-1-phosphate uridylyltransferase [OS=Saccha | 41% | [Progress] | 19 | 19 | 3.355 | 2.869 | 3.531 | 335.4 | 286.9 | 353.0 | 100.0 | 2.41 | 2.90 | 1.75 | 2.05 | [Green] | [Green] | [Green] | [Green] |
| ... | ... | ... | ... | ... | ... | ... | ... | ... | ... | ... | ... | ... | ... | ... | ... | ... | ... | ... | ... | ... |
| P31539 | heat shock protein 104 [OS=Saccharomyces cerevisiae S28c | 38% | [Progress] | 32 | 32 | 1.586 | 0.997 | 0.799 | 158.6 | 99.7 | 79.9 | 100.0 | 1.57 | 1.96 | 2.39 | 0.54 | [Green] | [Green] | [Green] | [Green] |
| P00427 | Cytochrome c oxidase subunit 6, mitochondrial [OS=Saccha | 13% | [Progress] | 1 | 1 | 1.567 | 1.501 | 1.260 | 156.7 | 150.1 | 126.0 | 100.0 | 2.98 | 0.79 | 2.80 | 3.26 | [Green] | [Green] | [Green] | [Green] |

Show Associated Tables

Ready 1012/1289 Proteins; 1012 Protein Groups; 5777 Peptide Groups; 6372 PSMs; 17246 MS/MS Spectrum Info; 17252 Quan Spectra; 1/2 Input Files; 2 Specialized Traces; 291 Result Statistics

Result Statistics




$\Delta ura2$ Over-represented Pathways:

- Reactome
 - Pyrimidine biosynthesis
 - Metabolism of amino acids and derivatives
 - Nucleobase biosynthesis
- KEGG
 - Biosynthesis of antibiotics
 - Alanine, aspartate and glutamate metabolism
 - Pyrimidine metabolism
 - Glycine, serine and threonine metabolism
 - Biosynthesis of secondary metabolites
 - One carbon pool by folate
 - Glyoxylate and dicarboxylate metabolism
- BioCyc
 - UMP biosynthesis II
 - urea cycle
 - Nitrogen Compounds Metabolism
 - Pyrimidine Ribonucleotides De Novo Biosynthesis
 - de novo biosynthesis of pyrimidine ribonucleotides

Comparison: $\Delta ura2$ and $\Delta met6$

| List name | Proteins |
|--------------------|----------|
| ● Met6ko ms3 Lumos | 33 |
| ● Ura2ko ms3 Lumos | 21 |




Lists Comparison Proteins in Set

«Unique to Ura2ko ms3 Lumos | Taxonomy: Saccharomyces cerevisiae S288C


Proteins in Set

| Acc. Key | Description | Gene(s) | Length |
|----------|---|-------------|--------|
| P0C2J0-1 | Transposon Ty1-PR2 Gag-Pol polyprotein | YPR158... 3 | 1756 |
| P07259 | Protein URA2 | URA2 3 | 2214 |
| P24031 | Constitutive acid phosphatase | PHO3 3 | 467 |
| P17064 | Purine-cytosine permease FCY2 | FCY2 3 | 533 |
| P28272 | Dihydroorotate dehydrogenase (Fumarate) | URA1 3 | 314 |
| P07273 | transcription elongation factor S-II | DST1 3 | 309 |
| P50861 | 6,7-dimethyl-8-ribityllumazine synthase | RIB4 3 | 169 |
| P20051 | dihydroorotase | URA4 3 | 364 |
| P40054 | D-3-phosphoglycerate dehydrogenase 1 | SER3 3 | 469 |


Detailed Result for Selected Protein- KO Protein Met6




Detailed Result for Selected Peptide Group from KO Protein Met6




Quan Results- Details peptides




Mass Accuracy



MS2 Injection Time, Front Optics, LC




Reporter Ions S/N; Quad Status




Example Plots for Lumos SPS, 50 min QC run

QC Tools Instrument: Housekeeping Proteins Ratios Should be a Perfect 1




- Scaled ,
- No normalization
- Correction factors applied


Mass Accuracy



MS2 Injection Time, Front Optics, LC



Reporter Ions S/N; Quad Status



Fusion SPS data, instrument needs maintenance

QC Tools Instrument: Housekeeping Proteins Ratios Should be a Perfect 1



Thermo Proteome Discoverer 2.2.0.388


File View Administration Tools Window Help

Study: TNO standard x TKOTT11_1ms3_2unique75 x TKOTT11_1ms3_2unique75(1) x Administration x TKOTT11_1ms3_2 x TKOTT11_1ms2_2_42 x 171003_TMT11_50min_HFX_SN6042 x TKOtm11_sps-1_1hral1 x TKOtm11_ms2al_1hral2(1) x

Proteins Protein Groups Peptide Groups PSMs MS/MS Spectrum Info Quan Spectra Input Files Specialized Traces Result Statistics

| | Checked | Master | Accession | Description | Coverage [%] | Sequence Coverage | # Peptides | # PSMs | # Unique Peptides | # Protein Groups | # AAs | MW [kDa] | calc. pI | Score S | # Peptides | Biological Process | Cellular Compo |
|----|-------------------------------------|-------------------------------------|-----------|---|--------------|-------------------|------------|--------|-------------------|------------------|-------|----------|----------|---------|------------|--------------------|----------------|
| 1 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P00925 | Enolase 2 [OS=Saccharomyces cerevisiae S288C] | 84% | | 28 | 64 | 9 | 1 | 437 | 46.9 | 6.00 | 255.28 | 28 | | |
| 2 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P00924 | Enolase 1 [OS=Saccharomyces cerevisiae S288C] | 84% | | 30 | 58 | 11 | 1 | 437 | 46.8 | 6.62 | 219.93 | 30 | | |
| 3 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P07149 | Fatty acid synthase subunit beta [OS=Saccharomyces cerevisiae S288C] | 36% | | 56 | 57 | 56 | 1 | 2051 | 228.5 | 5.92 | 212.20 | 56 | | |
| 4 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P19097 | Fatty acid synthase subunit alpha [OS=Saccharomyces cerevisiae S288C] | 37% | | 56 | 57 | 56 | 1 | 1887 | 206.8 | 5.44 | 211.39 | 56 | | |
| 5 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P32324 | Elongation factor 2 [OS=Saccharomyces cerevisiae S288C] | 67% | | 43 | 49 | 43 | 1 | 842 | 93.2 | 6.32 | 207.19 | 43 | | |
| 6 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P00549 | Pyruvate kinase 1 [OS=Saccharomyces cerevisiae S288C] | 71% | | 30 | 46 | 30 | 1 | 500 | 54.5 | 7.68 | 195.68 | 30 | | |
| 7 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P00359 | Glyceraldehyde-3-phosphate dehydrogenase 3 [OS=Saccharomyces cerevisiae S288C] | 86% | | 22 | 59 | 5 | 1 | 332 | 35.7 | 6.96 | 194.26 | 22 | | |
| 8 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P17255 | V-type proton ATPase catalytic subunit A [OS=Saccharomyces cerevisiae S288C] | 56% | | 50 | 50 | 50 | 1 | 1071 | 118.6 | 6.16 | 187.77 | 50 | | |
| 9 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P11484 | Ribosome-associated molecular chaperone SSB1 [OS=Saccharomyces cerevisiae S288C] | 72% | | 34 | 39 | 3 | 1 | 613 | 66.6 | 5.44 | 164.06 | 34 | | |
| 10 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P15108 | ATP-dependent molecular chaperone HSC82 [OS=Saccharomyces cerevisiae S288C] | 62% | | 39 | 42 | 8 | 1 | 705 | 80.8 | 4.83 | 163.99 | 39 | | |
| 11 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P10592 | Heat shock protein SSA2 [OS=Saccharomyces cerevisiae S288C] | 68% | | 34 | 40 | 10 | 1 | 639 | 69.4 | 5.06 | 161.11 | 34 | | |
| 12 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P40150 | Ribosome-associated molecular chaperone SSB2 [OS=Saccharomyces cerevisiae S288C] | 68% | | 33 | 38 | 2 | 1 | 613 | 66.6 | 5.49 | 159.25 | 33 | | |
| 13 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P16861 | ATP-dependent 6-phosphofructokinase subunit alpha [OS=Saccharomyces cerevisiae S288C] | 55% | | 40 | 40 | 39 | 1 | 987 | 107.9 | 6.39 | 156.38 | 40 | | |
| 14 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P00358 | glyceraldehyde-3-phosphate dehydrogenase 2 [OS=Saccharomyces cerevisiae S288C] | 82% | | 21 | 45 | 4 | 1 | 332 | 35.8 | 6.96 | 156.11 | 21 | | |
| 15 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P00560 | Phosphoglycerate kinase [OS=Saccharomyces cerevisiae S288C] | 77% | | 26 | 40 | 26 | 1 | 416 | 44.7 | 7.61 | 154.68 | 26 | | |
| 16 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P02829 | ATP-dependent molecular chaperone HSP82 [OS=Saccharomyces cerevisiae S288C] | 59% | | 37 | 40 | 6 | 1 | 709 | 81.4 | 4.91 | 154.09 | 37 | | |
| 17 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P06169 | pyruvate decarboxylase isozyme 1 [OS=Saccharomyces cerevisiae S288C] | 58% | | 22 | 35 | 20 | 1 | 563 | 61.5 | 6.19 | 148.89 | 22 | | |
| 18 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | P10591 | heat shock protein SSA1 [OS=Saccharomyces cerevisiae S288C] | 67% | | 33 | 37 | 9 | 1 | 642 | 69.6 | 5.11 | 145.00 | 33 | | |


Show Associated Tables




- Scaled ,
- No normalization
- Correction factors applied

Yeast TKO TMT 11 Standard: QC run Expected Average Instrument Performance


Protein Groups, 50 min run



Unique Peptides, 50 min run




Accuracy of Quantitation, 50 min runs




MF= SPS Mass Matches (%), 65%
New feature in PD 2.3

Yeast TKO TMT 11 Standard: QC run Expected Average Instrument Performance


Protein Groups, 50 min run



Unique Peptides, 50 min run




Accuracy of Quantitation, Lumos SPS 50 min run




Yeast TKO TMT 11 Standard 2 hour run – Method Development/Optimization

RT: 0.00 - 185.00





NL:
7.09E9
Base Peak
MS
TKOTT11_1ms2_2hr1


Protein Groups, 120 min gradient



P15992, Heat Shock protein 26(301ppm)




Unique Peptides, 120 min gradient




Summary: Product Information- Available October 1

- Excellent tool for LC and mass spec method development
- Excellent QC assay tool for quality assessment of LC and mass spec instrument status
- Provides accuracy, precision and dynamic range assessments for different mass spec strategies
- Optional: TMT labeled Peptide Retention Time Calibration mix (PRTC, yeast heavy isotope peptides) can be spiked in for triggered, multiplexed assay (TOMAHAQ) method development
- **A40938** Pierce TMT11plex yeast digest standard, 20µg
- **A40939** Pierce TMT11plex yeast digest standard, 5 x 20µg



Pierce TMT11 yeast TKO standard



A Standard Multiplexed Targeted Proteomic Assay Utilizing Isobaric Labels for Evaluation of TOMAHAQ

Christopher M. Rose¹, Rosa Viner², Jae Cho³, John Rogers⁴, Devin K. Schweppe⁴, Brian K. Erickson⁴, Steven P. Gygi⁴, Donald S. Kirkpatrick¹

¹Genentech Inc., South San Francisco, CA; ²Thermo Fisher Scientific, San Jose, CA; ³Thermo Fisher Scientific, Rockford, IL; ⁴Harvard Medical School, Boston, MA



Introduction

Triggered by Offet, Multiplexed, Accurate mass, High resolution, and Absolute Quantification (TOMAHAQ) is a recently introduced targeted proteomics method that combines peptide and sample multiplexing. Establishing TOMAHAQ methods within a lab has proven challenging due to the various formulations of TOMAHAQ experiments and lack of a common assay that can be universally implemented.

TOMAHAQ initially combined samples multiplexed with TMT10 and TMT0 labeled synthetic trigger peptides. However, as more labs establish TOMAHAQ assays using available reagents, many questions have been raised:

- 1) Should trigger peptides be light, heavy, or another protease;
- 2) How many peptides should be targeted for method development?
- 3) Is method failure due to synthetic peptides, cellular lysate, or instrument performance?

At the core of these questions is the lack of a standard TOMAHAQ assay that can be used to develop and establish TOMAHAQ methods within a laboratory. Here, we introduce a standard TOMAHAQ assay that utilizes commercially available heavy peptides and whole cell lysate that can be used to create a standard assay for the establishment and evaluation of TOMAHAQ within research laboratories. We combine this assay with recent updates to a freely available program, TomahqCompanion, which can be used to make TOMAHAQ instrument methods and analyze TOMAHAQ data.

TOMAHAQ + TomahqCompanion Workflow





Figure 2. TOMAHAQ sample analysis workflow using TomahqCompanion. There are two stages of analysis. (1) Printing run. This is an analysis of five peptides alone that is required when developing an assay. A target list for the printing run is generated by TomahqCompanion and the results of the printing run can be viewed by TomahqCompanion. The printing run can be combined with a separate TOMAHAQ method to produce a TOMAHAQ method file. (2) TOMAHAQ run. After generating the TOMAHAQ method you can run the synthetic peptides with your sample of interest and analyze the experimental sample. The final data can be viewed, sorted, and exported using TomahqCompanion.

A standard mixture for TOMAHAQ evaluation

















Figure 3. A standard mixture for TOMAHAQ evaluation. The standard mixture comprises Pierce RT calibration (PRTC) peptides, HeLa digest, TMT10 plus, and S10-TMT reagents. The PRTC peptides were labeled with S10-TMT and spiked into a mixture of HeLa digest (labeled with 8 channels of standard 1:1) and PRTC peptides labeled with 8 channels of TMT10 plus mixed 10 : 5 : 2.5 : 1 : 0 : 0 : 1 : 2.5 : 5 : 10. Characteristic TOMAHAQ analysis in a complex background. TMT10 plus labeled PRTC peptides were mixed at varying amounts relative to the TMT10 plus HeLa digest background.

Dilution series of standard mixture



Number of notches



TOMAHAQ performed through Instrument API

The current implementation of TOMAHAQ through the method editor allows ~150 peptides to be analyzed, but increasing the number of peptides a more advanced TOMAHAQ implementation needed to be created. The instrument API allows a program written in C# to receive data in real time and send scan instructions. Performing TOMAHAQ through the API will enable:

1. Increase the number of peptides that can be monitored into the thousands
2. Enable more advanced retention time scheduling (elution order scheduling)
3. Allow custom filters to be applied to increase the quantitative accuracy
4. Integration with TomahqCompanion program (also written in C#)



TOMAHAQ API Scan Sequence




Figure 19. TOMAHAQ performed through the Instrument API was used to analyze the TOMAHAQ standard run. Each of the peptides was identified and targeted through API directed analysis. The data was visualized using TomahqCompanion, which reveals that the API implementation is correctly targeting the complex TOMAHAQ scan structure.

Conclusions and future directions

- Developed a TOMAHAQ standard using components that can be either purchased (HeLa, Pierce RT calibration peptides) or acquired (TMT10) from Thermo
- Quantifying the best peptide (lowest amount of interference) demonstrates superior quantitative accuracy throughout all conditions tested
- Dilution of TOMAHAQ standard demonstrates quantitative accuracy maintained down to 100 amol in lowest channel (RE + 20% for best peptide)
- Decreasing the number of notches increases quantitative accuracy, but limits reporter ion signal in the MS3
- Longer runs increase quantitative accuracy, with the longest run tested (90 min) returning the best quantitative results
- Precursor range exclusion setting does not appear to have a large impact on quantitative accuracy
- TOMAHAQ implemented through the API successfully for TOMAHAQ standard. Future work will increase the number of peptides that can be quantified using more advanced RT scheduling and quantitative filters

Acknowledgments

We would like to thank the members of the Microchemistry, Proteomics, and Lipidomics department in Genentech Research & Early Development. We would also like to thank the Thermo research and development team, specifically Denis Bailey and Philip Barnes for their work on increasing the number of peptides that can be analyzed through the standard method editor.

References

Bailey, D.E., Barnes, P., Rose, C.M., Viner, R., Cho, J., Rogers, J., Schweppe, D.K., Erickson, B.K., Gygi, S.P., Kirkpatrick, D.S. (2017) *Journal of Proteomics*, 141, 1-12.

Acknowledgments

- Jae Choi
- Aaron M. Robitaille
- John C. Rogers
- Kay Opperman
- David Horn
- Devin Drew
- Tabi N. Arrey
- Andreas Huhmer