




ThermoFisher
S C I E N T I F I C

Intact and Top-Down Protein Analysis Guide: Orbitrap Fusion Lumos


Seema Sharma
October 2016



The world leader in serving science



Orbitrap Fusion Lumos





Application Dependency: Pressure & Resolution

MS
Protein
>50 kDa

MS
Protein
<40kDa

MS/MS


MS
Protein
40-50 kDa

Standard Pressure Mode
R @ 15000


Intact Protein Mode
IRM @ 2mtorr
R @ 120000

Intact Protein Mode
IRM @ 2mtorr
R @ 120000


Intact Protein Mode
IRM @ 1mtorr
R @ 240000




Intact mAb Charge States



Carbonic Anhydrase Isotopic Resolution




Carbonic Anhydrase HCD MS/MS





Reduced mAb Heavy Chain Isotopic Resolution

General Calibrations



- Check all positive calibrations before proceeding to extended positive calibrations.
- Check Negative mode Ion optics & Ion Trap electron multiplier gain calibrations prior to ETD calibration.



Extended Calibrations

MS
Protein
>50 kDa

Normal Mass Range
No Extended
Calibrations Required

High Mass Range
Positive Extended:
High Mass Range
CalMix & Enfuvirtide

- ▶ High Mass Range : Calibration Solution
- ▶ High Mass Range : Enfuvirtide

Extended Calibrations Contd.

MS
Protein
<40 kDa

Set IRM Pressure @
2mtorr
Positive Extended:
Intact Protein Mode

Intact Protein Mode : Calibration Solution

Optics

- Source to IRM
- IRM to C-Trap
- IRM to Ion Trap

MS
Protein
40-50 kDa

Set IRM Pressure @
1mtorr
Positive Extended:
Intact Protein Mode

Intact Protein Mode : Calibration Solution

Optics

- Source to IRM
- IRM to C-Trap
- IRM to Ion Trap

IPM Pressure Setting

- IPM pressure can be changed in diagnostics, by changing the parameter value under 'Set IRM pressure (Intact)' and then click on Start.
- Hardware → System → IRM → Set IRM Pressure (Intact) → 'IRM_Pressure' Parameter Value
- Note: The instrument can only be calibrated for one 'IRM Pressure (Intact)' at a given time.
- The intact protein mode calibration (optics) needs to be run at the desired IRM pressure (Intact).

The screenshot displays the 'Diagnostics' menu in the software. The 'System' folder is expanded, showing 'IRM' selected. Under 'IRM', 'Set IRM Pressure (Intact)' is highlighted. Below this, there are checkboxes for 'Check IRM Pressure Regulation' and 'Adjust Collision Pressure Readback'. The 'IRM_Pressure' parameter is set to 0.001. At the bottom right, a 'Start' button is visible.

Parameter Name	Parameter Value
IRM_Pressure	0.001

Extended Calibrations Contd.

MS/MS

Set IRM Pressure @
2mtorr

Positive Extended:
Intact Protein Mode
&
ETD Extended:
Intact Protein Mode

Intact Protein Mode : Calibration Solution

Optics

- Source to IRM
- IRM to C-Trap
- IRM to Ion Trap


ETD Extended

Intact Protein Mode

- Source to IRM

FT Signal: Averaged Transients vs. Averaged Scans

- To improve spectra S/N use either microscans (each is a single transient) or Averager in Tune. S/N improves as square root of number of transients averaged.
- When spectra are averaged in Qual browser, both signal and noise are averaged, the gain in S/N is somewhat lower.
- Whenever possible try to average as many transients as possible to improve spectral S/N.
- Averager updates the scan with every transient added.
- Microscans update the scan only when all transients are added.
- The end result is exactly the same whether you use n microscans or average n transients through Averager. The only difference is in the scan update frequency.
- Averager uses rolling average routine.
- When using microscans (n) the scan appear after averaging the last n transients.



ThermoFisher
S C I E N T I F I C

Direct Infusion: Carbonic Anhydrase


The world leader in serving science

Carbonic Anhydrase (29 kDa)

- P/N C2522, Sigma
- Direct Infusion: Flow rate 3ul/min
- 1-5 pmol/ul in 50:50 AcN: H₂O, 0.1% FA
- HESI probe (same source conditions as for calmix)

(42.0106) SHHWGYGKH NGPEHWHKDF PIANGERQSP
VDIDTKAVVQ DPALKPLALV YGEATSRRMV NNGHSFNVEY
DDSQDKAVLK DGPLTGTYRL VQFHFHWGSS DDQGSEHTVD
RKKYAAELHLVHWNTKYGDF GTAAQQPDGL AVVGVFLKVG
DANPALQKVL DALDSIKTKG KSTDFPNFDPGSLLPNVLDY
WTYPGSLTTP PLLESVTWIV LKEPISVSSQ QMLKFRTLNF
NAEGPELLM LANWRPAQPL KNRQVRGFPK

N-terminal
acetylation







Carbonic Anhydrase Full MS

240K, IPM 1mtorr

AGC 2e5, 20uscans




CA_1uscan_240K_1_mtorr_2e5_RF30_XT_00001_M_#2 RT: 2.00 AV: 1 NL: 3.59E7
T: FTMS + p ESI Full ms [400.00-1600.00]



Carbonic Anhydrase MS/MS CID

120K, 2mtorr, CID CE 25

AGC 5e5, 200 uscans



ProSight Lite

File Modify Experimental Data Modify Candidate Sequence Export Help About

N S H H W G Y G K H N G P E H W H K D F P I A N G E 25
26 R Q S P V D I D T T K A V V Q D P A L K P L A L V Y 50
51 G E A T S R R M V N N G H S F N V E Y D D S Q D K 75
76 A V L K D G P L T G T Y R L V Q F H F H W G S S D 100
101 D Q G S E H T V D R K K Y A A E L H L V H W N T K 125
126 Y G D F G T A A Q Q P D G L A V V G V F L K V G D 150
151 A N P A L Q K V L D A L D S I K T K G K S T D F P 175
176 N F D P G S L L P N V L D Y W T Y P I G S L L T T P P 200
201 L L L E S V T W I V L K E P I S V S S Q Q M L K F R 225
226 T L N F N A E G E P E L L L M L A N W R P A Q L P L K 250
251 N R Q V R G F P K C

Precursor Mass
Type: Monoisotopic
Observed: 29,006.68
Theoretical: 29,006.68
Mass Diff. (Da): 0.00
Mass Diff. (ppm): 0.02

Scores
PCS: 593.55
P-Score: 1.2e-53
% Fragments Expl... 24 %
% Residue Cleava... 33 %


Modification (S1)
No Modification
Custom 42.0106
Common

ProSight Lite:
10 ppm fragment mass
tolerance

Carbonic Anhydrase MS/MS HCD

120K, 2mtorr, HCD CE 10

AGC 5e5, 200 uscans



ProSight Lite

File Modify Experimental Data Modify Candidate Sequence Export Help About

N **S**H|H|W|G|Y G K|H|N G P E H|W|H|K|D|F|P|I|A|N|G|E 25
26 |R|Q|S|P|V|D|I|D|T|K A|V|V|Q D|P|A|L|K P L A|L|V Y 50
51 G E A T S R R M V N N G H S|F|N|V|E|Y D D|S Q D K 75
76 A V L K D G P L T G T Y R L V Q F H F H|W G S S D 100
101 D Q G S E H T V D R K K Y A A E|L H L V H W N T K 125
126 Y G D F G T|A|A|Q|Q|P D G|L A V V G V F L K V G D 150
151 A N P A L Q K V L D A L D S I K T K G K S T D F P 175
176 N F D P G S L L|P N V L D Y|W T|Y|P G S|L|T|T|P|P 200
201 L L L E S V|T|W|I|V L|K|E|P I|S V S S Q Q M L K F R 225
226 T L N F N|A E G E|P|E|L|L|L|L|L|A|N|W|R|P A Q|P L K 250
251 N R Q V R G F P K C

Precursor Mass
Type: Monoisotopic
Observed: 29,006.68
Theoretical: 29,006.68
Mass Diff. (Da): 0.000
Mass Diff. (ppm): 0.02

Scores
PCS: 609.77
P-Score: 6.4e-55
% Fragments Expl... 17 %
% Residue Cleava... 33 %


Modification (S1)
 No Modification
 Custom 42.0106
Common



Carbonic Anhydrase ETD HD

120K, 2mtorr, ETD Rxn 10ms

AGC 1e6, Reag AGC 7e5, 200 uscans



ProSight Lite

File | Modify Experimental Data | Modify Candidate Sequence | Export | Help About

N **S** H H W G Y G K H N G P E H W H K D F P I A N G E 25
 26 R Q S P V D I D T K A V V Q D P A L K P L A L V Y 50
 51 G E A T S R R M V N N G H S F N V E Y D D S Q D K 75
 76 A V L K D G G P L T G T Y R L V Q F H F H W G S S D 100
 101 D Q G S E H T V D R K K Y A A E L H L V H W N T K 125
 126 Y G D F G T A A Q Q P D G L L A V V G L V F L K L V G 150
 151 A N P A L Q K V L D A L D S I K T K G K S L T D F P 175
 176 N F D P G S L L P N V L L D Y W T Y P G S L T T P P 200
 201 L L E S V T W I V L K E P I S V S S Q Q M L K L F R 225
 226 T L N F N A E G E P E L L M L L A N W R P A Q P L K 250
 251 N R Q V R G F P K C

Precursor Mass
 Type: Monoisotopic
 Observed: 29,006.68
 Theoretical: 29,006.68
 Mass Diff. (Da): 0.000
 Mass Diff. (ppm): 0.02

Scores
 PCS: 1,132.10
 P-Score: 5.5e-95
 % Fragments Expl: 31.8%
 % Residue Cleava... 51%

Modification (S1)
 No Modification
 Custom 42.0106
 Common



ThermoFisher
SCIENTIFIC

Direct Infusion: Enolase



The world leader in serving science

Enolase (47 kDa)

- P/N E6126, Sigma
- Direct Infusion: Flow rate 3ul/min
- Clean-up using Bio-Rad Micro spin columns (follow manufacturer's instructions)
- 10 -20 pmol/ul in 50:50 AcN: H₂O, 0.1% FA

AVSKVYARSVYDSRGNPTVEVELTTEKGVFRSIVPSGASTGVHEALEMRDGDKSKWMGKG
VLHAVKNVNDVIAPAFVKANIDVKDQKAVDDFLISLDGTANKSKLGANAILGVSLAASRAAAAE
KNVPLYKHLADLSKSKTSPYVLPVPFLNVLNGGSHAGGALALQEFMIAPTGAKTFAEALRIGS
EVYHNLKSLTKKRYGASAGNVGDEGGVAPNIQTAEALDLIVDAIKAAGHDGKVKIGLDCASS
EFFKDGKYDLDFKNPNSDKSKWLTGPQLADLYHSLMKRYPIVSIEDPFAEDDWEAWSHFFKT
AGIQIVADDLTVTNPKRIATAIEKKAADALLKVNQIGTLSESIKAAQDSFAAGWGVMVSHRSG
ETEDTFIADLVVGLRTGQIKTGAPARSERLAKLNQLLRIEEELGDNAVFAGENFHHGDKL

Enolase Full MS







Enolase 47 kDa, MS/MS CID

240K, 1mtorr, CID CE 35


AGC 5e5, 200 uscans







Enolase MS/MS HCD




Enolase ETD HD CS 52+







Enolase CS 52+, ETD HD Sequence Coverage



Enolase ETD HD CS 63+



Enolase CS 63+, ETD HD Seq Coverage




Enolase Combined Sequence Coverage

56% Seq Coverage

```
N [A] V S K[V] Y[A] R[S] V[Y] D[S] R[G] N[P] T[V] E[V] E[L] T[T] 25
26 [E] K[G] V[F] R[S] I[V] P[S] G A[S] T[G] V H[E] A[L] E[M] R[D] 50
51 [G] D[K] S[K] W[M] G[K] G[V] L[H] A V K[N] V[N] D V[I] A[P] A 75
76 [F] V[K] A[N] I[D] V[K] D[Q] K[A] V[D] D[F] L[I] S[L] D G[T] A 100
101 [N] K S K[L] G[A] N[A] I[L] G[V] S[L] A[A] S R A A A[A] E K 125
126 N V P L[Y] K H L A[D] L S[K] S[K] T S P Y V[L] P[V] P[F] 150
151 [L] N[V] L[N] G G S H A G G[A] L[A] L[Q] E[F] M[I] A P T G 175
176 A K T F A E A L R I G S E V Y H[N] L K S L T K K R 200
201 Y G A S A G N V[G] D[E] G[G] V[A] P N I Q[T] A[E] E A L 225
226 [D] L I V D A I K A A G H D G K I K I G L D C A S S 250
251 [E] F F[K] D G K Y D L D F K[N] P[N] S D K[S] K W L T G 275
276 P Q L A D L Y H S L M K R[Y] P I V S I E D[P] F[L] A[E] 300
301 [D] D W E A W S H[F] F[K] T A G I Q I[V] A[D] D L[T] V T 325
326 [N] P K[R] I A T[A] I[E] K[K] A[A] D[A] L L L[K] V[N] Q I G 350
351 [T] L[S] E[S] I K[A] A Q[D] S[F] A A[G] W[G] V[M] V S H[R] S 375
376 G[E] T[E] D[T] F[I] A[D] L[V] V[G] L R[T] G[Q] I[K] T G[A] P 400
401 [A] R[S] E[R] L[A] K[L] N[Q] L[L] R[I] E[E] E[L] G[D] N[A] V[F] 425
426 [A] G[E] N[F] H[H] G D K L C
```

Enolase MS2 CID MS3 CID



Enolase MS2 CID MS3 ETD

MS2 CID m/z 898, MS3 ETD m/z 870

ProSight Lite

File Modify Experimental Data Modify Candidate Sequence Export Help About


N **A** V S K[V]Y[A]R[S]V[Y]D[S]R[G]N P[T]V[E]V[E]L[T]T 25
26[E]K[G]V F[R]S[I]V P[S]G A[S]T[G]V H[E]A L[E]M[R]D 50
51[G]D[K]S K[W]M G[K]G V[L]H A[V]K[L]N[V]N[D]V I A P A 75
76 F V K A[N]I D V K[D]Q K A V D[D]F L I[S] L D G T A 100
101 N K[S]K L G A I N A I L G V S L A A S R A A A I A E K 125
126 N V P L Y K H L A D L S K[S]K T S P Y V L P V C

MS3 Precursor Sequence

Precursor Mass
Type: Monoisotopic
Observed: 46,642.21
Theoretical: 15,637.31
Mass Diff. (Da): 31,004.905
Mass Diff. (ppm): n/a

Scores
PCS: 605.81
P-Score: 1.3e-54
% Fragments Expl... 37 %
% Residue Cleava... 36 %

Modification (V148)
No Modification
Custom -18.01





ThermoFisher
SCIENTIFIC

DDA LC-MS Analysis: Bacterial Lysates

The world leader in serving science

Chromatography

- LC: Dionex Ultimate 3000RSLCnano
- Flow rate: 12 $\mu\text{l}/\text{min}$
- Column: ProSwift Monolithic RP-4H
200 μm x 25cm
- Trap Column (optional): Prototype C4, 200 μm x 5mm
- Solvents: A: H₂O, 0.1% FA
B: AcN, 0.1% FA
- Gradient: 20% B to 60% B in 90min




Instrument Method Full MS

MS Scan Properties	
Detector Type	Orbitrap
Orbitrap Resolution	120000
Mass Range	Normal
Use Quadrupole Isolation	<input checked="" type="checkbox"/>
Scan Range (m/z)	500-1700
RF Lens (%)	30
AGC Target	5.0e5
Maximum Injection Time (ms)	100
Microscans	5
Data Type	Profile
Polarity	Positive
Source Fragmentation	<input type="checkbox"/>
Use EASY-IC	<input type="checkbox"/>

- High / High Method
 - Full MS @ 120000 resolution
 - MS2 @ 120000 resolution
- Medium / High Method
 - Full MS @ 15000 resolution
 - MS2 @ 120000 resolution

Instrument Method DDA



Charge State Properties	
Include charge state(s)	5-24
Include undetermined charge states	<input checked="" type="checkbox"/>
Include charge states 25 and higher	<input checked="" type="checkbox"/>

Dynamic Exclusion Properties	
Exclude after n times	1
Exclusion duration (s)	60
Mass Tolerance	<input type="radio"/> ppm <input checked="" type="radio"/> m/z
Low	2.00
High	2.00

Decisions Properties	
Data dependent mode	<input type="radio"/> Top Speed <input checked="" type="radio"/> Top N
Precursor Priority	Most Intense
Number of Scan Event Types : 1	
SCAN EVENT TYPE : 1	
Intensity greater than	5.0e5


DDA Parameters

Data-Dependent MS ⁿ Scan Properties	
MS ⁿ Level	2
Isolation Mode	Quadrupole
Isolation Window (m/z)	5
Use Isolation m/z Offset	<input type="checkbox"/>
Activation Type	ETD
Use Calibrated Charge Dependent ETD Parameters	<input type="checkbox"/>
ETD Reaction Time (ms)	10
ETD Reagent Target	7.0e5
Maximum ETD Reagent Injection Time (ms)	200
ETD Supplemental Activation	<input type="checkbox"/>
Detector Type	Orbitrap
Scan Range Mode	Auto: m/z Normal
Orbitrap Resolution	120000
First Mass (m/z)	200
AGC Target	5.0e5
Inject Ions for All Available Parallelizable Time	<input checked="" type="checkbox"/>
Maximum Injection Time (ms)	400
Microscans	5
Data Type	Profile

Data-Dependent MS ⁿ Scan Properties	
MS ⁿ Level	2
Isolation Mode	Quadrupole
Isolation Window (m/z)	5
Use Isolation m/z Offset	<input type="checkbox"/>
Activation Type	CID
CID Collision Energy (%)	30
Activation Q	0.25
Multistage Activation	<input type="checkbox"/>
Detector Type	Orbitrap
Scan Range Mode	Auto: m/z Normal
Orbitrap Resolution	120000
AGC Target	5.0e5
Inject Ions for All Available Parallelizable Time	<input checked="" type="checkbox"/>
Maximum Injection Time (ms)	400
Microscans	5
Data Type	Profile
Use EASY-IC	<input type="checkbox"/>

Data-Dependent MS ⁿ Scan Properties	
MS ⁿ Level	2
Isolation Mode	Quadrupole
Isolation Window (m/z)	5
Use Isolation m/z Offset	<input type="checkbox"/>
Activation Type	HCD
HCD Collision Energy (%)	10
Stepped Collision Energy	<input type="checkbox"/>
Detector Type	Orbitrap
Scan Range Mode	Auto: m/z Normal
Orbitrap Resolution	120000
First Mass (m/z)	200
AGC Target	5.0e5
Inject Ions for All Available Parallelizable Time	<input checked="" type="checkbox"/>
Maximum Injection Time (ms)	400
Microscans	5
Data Type	Profile
Use EASY-IC	<input type="checkbox"/>

E. coli DDA Total Ion Chromatogram





Data Processing: ProSight PD

Proteins											Peptide Groups	PSMs	MS/MS Spectrum Info
Master	Accession	Description	# PSMs	# Protein Groups	# AAs	MW [kDa]	calc. pI	Entrez Gene I	Gene ID	Biological Process	Cellular Component	Molecular Function	
✓	P0AC92	Protein GnsA	2	1	57	6.6	5.40		ER3413_101	cell communication cell death cell differentiation cell division cell growth cell organization and biogenesis cell proliferation cellular component movement cellular homeostasis coagulation defense response development metabolic process regulation of biological process reproduction response to stimulus transport	cell surface chromosome cytoplasm cytoskeleton cytosol endoplasmic reticulum endosome extracellular Golgi membrane nucleus organelle lumen proteasome ribosome spliceosomal complex vacuole	antioxidant activity catalytic activity DNA binding enzyme regulator activity metal ion binding motor activity nucleotide binding protein binding receptor activity RNA binding signal transducer activity structural molecule activity transcription regulator activity translation regulator activity transporter activity	
✓	P0AAZ7	UPF0434 protein YcaR	6	1	60	6.9	5.01		ER3413_938				
✓	P69913	Carbon storage regulator	7	1	61	6.9	8.62		ER3413_277				
✓	P0AEG8	Protein DsrB	2	1	62	6.9	4.65		ER3413_200				
✓	P0A7M6	50S ribosomal protein L29	21	1	63	7.3	9.99		ER3413_339				
✓	P0A8H8	DNA gyrase inhibitor YacG (ECO:0000255 HAMAP-Rule:MF_01217)	2	1	65	7.3	4.56		ER3413_97				
✓	P0A9Y6	Cold shock-like protein CspC	23	1	69	7.4	7.24		ER3413_187				
✓	P0A9X9	Cold shock protein CspA	6	1	70	7.4	5.95		ER3413_365				
✓	P0A972	Cold shock-like protein CspE	12	1	69	7.5	8.54		ER3413_637				
✓	Q2M7R5	Uncharacterized protein YibT	10	1	69	8.0	9.41		ER3413_371				
✓	P0A8R4	Protein StyX	2	1	72	8.2	4.92		ER3413_343				
✓	P0AD07	Uncharacterized protein YecF	2	1	74	8.2	5.31		ER3413_196				
✓	P0AD24	UPF0352 protein YejL (ECO:0000255 HAMAP-Rule:MF_0000255)	8	1	75	8.3	5.78		ER3413_225				
✓	P68206	UPF0337 protein YjbJ	7	1	69	8.3	5.55		ER3413_417				
✓	P0AB14	Uncharacterized protein YccJ	18	1	75	8.5	4.74		ER3413_102				
✓	P0AA31	UPF0033 protein YedF	1	1	77	8.6	5.01		ER3413_198				
✓	P0A6A8	Acyl carrier protein (ECO:0000255 HAMAP-Rule:MF_01217)	6	1	78	8.6	4.06		ER3413_112				
✓	P0ACW6	Uncharacterized protein YdcH	3	1	74	8.9	9.31		ER3413_145				
✓	P0A8G9	Exodeoxyribonuclease 7 small subunit	4	1	80	8.9	4.41		ER3413_431				
✓	P0A7T7	30S ribosomal protein S18	1	1	75	9.0	10.59		ER3413_434				
✓	P0AA04	Phosphocarrier protein HPr	19	1	85	9.1	5.83		ER3413_249				
✓	P0AD10	Uncharacterized protein YecJ	7	1	83	9.1	5.05		ER3413_195				
✓	P0AC62	Glutaredoxin-3	3	1	83	9.1	7.24		ER3413_372				
✓	P0A7T3	30S ribosomal protein S16 (ECO:0000255 HAMAP-Rule:MF_01217)	14	1	82	9.2	10.55		ER3413_269				
✓	P0ACF4	DNA-binding protein HU-beta	25	1	90	9.2	9.70		ER3413_450				
✓	P0AB61	Protein YgiN	10	1	83	9.4	5.73		ER3413_129				
✓	P0A9W6	Acid stress protein IbaG (ECO:0000305)	4	1	84	9.4	6.28		ER3413_328				
✓	P0ACF0	DNA-binding protein HU-alpha	25	1	90	9.5	9.58		ER3413_412				
✓	P0AF36	Cell division protein ZapB	20	1	81	9.6	4.67		ER3413_404				
✓	P0AG63	30S ribosomal protein S17 (ECO:0000255 HAMAP-Rule:MF_01217)	1	1	84	9.7	9.60		ER3413_339				
✓	P0A8J4	UPF0250 protein YbeD	5	1	87	9.8	5.76		ER3413_645				
✓	P0AAN9	Anti-adaptor protein IraP	10	1	86	9.9	4.92		ER3413_390				
✓	P64540	Uncharacterized protein YfcL	1	1	92	10.0	4.41		ER3413_240				
✓	P0A800	DNA-directed RNA polymerase subunit omega	12	1	91	10.2	4.92		ER3413_376				

Protein Identifications



Higher number of proteins identified using High/High method





Identification of Large Proteins


Greater number of large proteins (>30kDa) identified using Med/High






High Flow LC-MS Analysis: Intact mAb

Intact IgG: Seven Major Glycosylated Forms




Deconvoluted Mass, Da	Reported Mass, Da
148018.5	
148220.6	148220.4
148383.0	148382.5
148546.0	148544.6
148709.6	148706.7
148871.2	148868.8
149035.7	



Why mass shifts?


How Many Components Are In Each Peak?



49+


- Covalent mods can include: SS red/ox, oxidation, deamidation, aa substitutions, Lys variants, sialated forms
- Non-covalent mods can include: metal adducts, phosphate adducts

Commonly Observed CS Distribution.....s????




Evidence For Multiple Overlapping Distributions

T: FTMS + p NSI sid=40.00 Full ms [1000.0000-5000.0000]




I: FTMS + p NSI sid=40.00 Full ms [1000.0000-5000.0000]



41+: Higher Resolution Reveals Multiple Isoforms

Unique to Orbitrap Analyzer



Chromatography

- Sigma mAb (p/n MSQC4)
- Sample Load: 1ug
- LC : Dionex Ultimate 3000
- Column: MabPac RP 2.1 x 10 cm

LC Gradient

Retention [min]	Flow [ml/min]	%B
0.000	0.300	10.0
0.000	0.300	10.0
2.000	0.300	30.0
7.000	0.300	45.0
8.000	0.300	90.0
9.000	0.300	10.0
10.000	0.300	10.0



Instrument Method Full MS

MS 01

Ion Source Properties

Ion Source Type	H-ESI ▾
Spray Voltage	<input checked="" type="radio"/> Static <input type="radio"/> Time Dependent
Positive Ion (V)	3500.00
Negative Ion (V)	3000.00
Sheath Gas (Arb)	60
Aux Gas (Arb)	20
Sweep Gas (Arb)	0
Ion Transfer Tube Temp (°C)	400
Vaporizer Temp (°C)	200
APPI Lamp	Not in Use ▾

MS Scan Properties

Detector Type	Orbitrap ▾
Orbitrap Resolution	15000 ▾
Mass Range	High ▾
Use Quadrupole Isolation	<input type="checkbox"/>
Scan Range (m/z)	2000-3500
RF Lens (%)	60
AGC Target	2.0e5 ▾
Maximum Injection Time (ms)	100
Microscans	10
Data Type	Profile ▾
Polarity	Positive ▾
Source Fragmentation	<input checked="" type="checkbox"/>
Energy (V)	50
Use EASY-IC	<input type="checkbox"/>

Divert Valve Set-up

Global Parameters

1.3 2.7 4 5.3 6.7 8

Ion Source

MS


Syringe

Divert Valve A


Divert Valve A Properties

	Time (min)	Position
1	0	1-2
2	1.5	1-6

Sigma mAb LCMS: 15K Full MS



BioPharma Finder Deconvolution Results



Average Mass	Matched Delta Mass (ppm)
146821.99	4.9
146660.35	8.4
146984.05	4.4
147146.78	8.4
146511.68	9.2



ThermoFisher
SCIENTIFIC

Targeted LC-MS Analysis: Reduced mAb

Chromatography: Flow Rate 300 ul/min

- Sample Load: 3ug
- Sigma mAb (p/n MSQC4)
Reduced deglycosylated
- LC : Dionex Ultimate 3000
(High Flow)
- Column: MabPac RP 2.1 x 10 cm

LC Gradient

Retention [min]	Flow [ml/min]	%B
0.000	0.300	10.0
0.000	0.300	10.0
2.000	0.300	30.0
7.000	0.300	45.0
8.000	0.300	90.0
9.000	0.300	10.0
10.000	0.300	10.0



Divert Valve Set-up

Global Parameters

1.3 2.7 4 5.3 6.7 8

3.5 (kV)
0 (kV)
-3.5 (kV)

Ion Source

MS

Syringe

Divert Valve A

1-2 1-6



Divert Valve A Properties

Import Export + X


	Time (min)	Position	
1	0	1-2	
2	1.5	1-6	

000


Deglycosylated Reduced mAb



Heavy Chain: ReSpect Sliding Window Deconvolution




Light Chain Isotopic Resolution




Light Chain: Xtract Deconvolution

deglyco_120K_run02_XT_00001_M_#2 RT: 2.00 AV: 1 NL: 2.85E8
T: FTMS + p ESI Full ms [600.00-1700.00]



Light Chain: Wide Isolation Window




Light Chain: ETD & EThcD

65% Sequence Coverage


N	Q	S	A	L	T	Q	P	R	S	V	S	G	S	P	G	Q	S	V	T	I	S	C	T	G	T	25	
26	S	S	D	I	G	G	Y	N	F	V	S	W	Y	Q	Q	H	P	G	K	A	P	K	L	M	I	50	
51	Y	D	A	T	K	R	P	S	G	V	P	D	R	F	S	G	S	K	S	G	N	T	A	S	L	75	
76	T	I	S	G	L	Q	A	E	D	E	A	D	Y	Y	C	C	S	Y	A	G	D	Y	T	P	G	100	
101	V	V	F	G	G	G	T	K	L	T	V	L	G	Q	P	K	A	A	P	S	V	T	L	F	P	125	
126	P	S	S	E	E	L	Q	A	N	K	A	T	L	V	C	L	L	I	S	D	F	Y	P	G	A	V	150
151	T	V	A	W	K	A	D	S	S	P	V	K	A	G	V	E	T	T	T	P	S	K	Q	S	N	175	
176	N	K	Y	A	A	S	S	Y	L	S	L	T	P	E	Q	W	K	S	H	R	S	Y	S	C	Q	200	
201	V	T	H	E	G	S	T	V	E	K	T	V	A	P	T	E	C	S	C								

- EThcD: 15ms, 25ms
- ETD iw 400: 7ms, 15ms

Heavy Chain: ETD Fragmentation



Heavy Chain:CID, HCD Fragmentation



Combined Seq Cov: ETD, HCD, CID

25% Sequence Coverage


- ETD 10ms, 15ms, 25ms
- CID 35%, HCD 10%

```
N  E V Q L V E S G G G L V Q P G G S L R L S C V A S 25
26  G F T L N N Y D M H W V R Q G I G K G L E W V S K 50
51  I G T A G D R Y Y A G S V K G R F T I S R E N A K 75
76  D S L Y L Q M N S L R V G D A A V Y Y C A R G A G 100
101 R W A P L G A F D I W G Q G T M V T V S S A S T K 125
126 G P S V F P L A P S S K S T S G G T A A L G C L V 150
151 K D Y F P E P V T V S W N S G A L T S G V H T F P 175
176 A V L Q S S G L Y S L S S V V T V P S S S L G T Q 200
201 T Y I C N V N H K P S N T K V D K K V E P K S C D 225
226 K T H T C P P C P A P E L L G G P S V F L F P P K 250
251 P K D T L M I S R T P E V T C V V V D V S H E D P 275
276 E V K F N W Y V D G V E V H N A K T K P R E E Q Y 300
301 N S T Y R V V S V L T V L H Q D W L N G L K E Y K C 325
326 K V S N K A L P A P I E K T I S K A K G Q P R E P 350
351 Q V Y T L P P S R D E L L T K N Q V S L T C L V K G 375
376 F Y P S D I A V E W E S N G Q P E N N Y K T T P P 400
401 V L D S D G S F F L Y S K L T V D K S R W Q Q G L N 425
426 V F S C S V M H E A L L H N H Y T Q K S L S L S P G C
```

Heavy Chain Isotopic Resolution: LC-MS

240K, 1mtorr, AGC 2e5

8ug sample on column, 5uscans, spectrum averaged over LC-MS peak



Chromatography: Flow Rate 15ul/min

- Sample load: 500ng
- Sigma mAb (p/n MSQC4)
Reduced deglycosylated
- LC : Dionex Ultimate 3000
(Microflow)
- Column: ProSwift Monolithic RP-4H
200um x 25cm



LC Gradient

Retention [min]	Flow [μ l/min]	%B
0.000	15.000	5.0
0.000	15.000	5.0
5.000	15.000	10.0
7.000	15.000	15.0
17.000	15.000	45.0
20.000	15.000	65.0
21.000	15.000	90.0
23.000	15.000	5.0
40.000	15.000	5.0

Light Chain: Top-down Fragmentation

HCD 10: 22% Seq Cov

N	Q	S	A	L	T	Q	P	R	S	V	S	G	S	P	G	Q	S	V	T	I	S	C	T	G	T	25	
26	S	S	D	I	G	G	Y	N	F	V	S	W	Y	Q	Q	H	P	G	K	A	P	K	L	M	I	50	
51	Y	D	A	T	K	R	P	S	G	V	P	D	R	F	S	G	S	K	S	G	N	T	A	S	L	75	
76	T	I	S	G	L	Q	A	E	D	E	A	D	Y	Y	C	C	S	Y	A	G	D	Y	T	P	G	100	
101	V	V	F	G	G	G	T	K	L	T	V	L	G	Q	P	K	A	A	P	S	V	T	L	F	P	125	
126	P	S	S	E	E	L	Q	A	N	K	A	T	L	V	C	L	L	I	S	D	F	Y	P	G	A	V	150
151	T	V	A	W	K	A	D	S	S	P	V	K	A	G	V	E	T	T	T	P	S	K	Q	S	N	175	
176	N	K	Y	A	A	S	S	Y	L	S	L	L	T	P	E	Q	W	K	S	H	R	S	Y	S	C	Q	200
201	V	T	H	E	G	S	T	V	E	K	T	V	A	P	T	E	C	S	C								

CID 30: 22% Seq Cov

N	Q	S	A	L	T	Q	P	R	S	V	S	G	S	P	G	Q	S	V	T	I	S	C	T	G	T	25	
26	S	S	D	I	G	G	Y	N	F	V	S	W	Y	Q	Q	H	P	G	K	A	P	K	L	M	I	50	
51	Y	D	A	T	K	R	P	S	G	V	P	D	R	F	S	G	S	K	S	G	N	T	A	S	L	75	
76	T	I	S	G	L	Q	A	E	D	E	A	D	Y	Y	C	C	S	Y	A	G	D	Y	T	P	G	100	
101	V	V	F	G	G	G	T	K	L	T	V	L	G	Q	P	K	A	A	P	S	V	T	L	F	P	125	
126	P	S	S	E	E	L	Q	A	N	K	A	T	L	V	C	L	L	I	S	D	F	Y	P	G	A	V	150
151	T	V	A	W	K	A	D	S	S	P	V	K	A	G	V	E	T	T	T	P	S	K	Q	S	N	175	
176	N	K	Y	A	A	S	S	Y	L	S	L	L	T	P	E	Q	W	K	S	H	R	S	Y	S	C	Q	200
201	V	T	H	E	G	S	T	V	E	K	T	V	A	P	T	E	C	S	C								

ETD 15ms: 46% Seq Cov

N	Q	S	A	L	T	Q	P	R	S	V	S	G	S	P	G	Q	S	V	T	I	S	C	T	G	T	25	
26	S	S	D	I	G	G	Y	N	F	V	S	W	Y	Q	Q	H	P	G	K	A	P	K	L	M	I	50	
51	Y	D	A	T	K	R	P	S	G	V	P	D	R	F	S	G	S	K	S	G	N	T	A	S	L	75	
76	T	I	S	G	L	Q	A	E	D	E	A	D	Y	Y	C	C	S	Y	A	G	D	Y	T	P	G	100	
101	V	V	F	G	G	G	T	K	L	T	V	L	G	Q	P	K	A	A	P	S	V	T	L	F	P	125	
126	P	S	S	E	E	L	Q	A	N	K	A	T	L	V	C	L	L	I	S	D	F	Y	P	G	A	V	150
151	T	V	A	W	K	A	D	S	S	P	V	K	A	G	V	E	T	T	T	P	S	K	Q	S	N	175	
176	N	K	Y	A	A	S	S	Y	L	S	L	L	T	P	E	Q	W	K	S	H	R	S	Y	S	C	Q	200
201	V	T	H	E	G	S	T	V	E	K	T	V	A	P	T	E	C	S	C								

ETD 25ms: 44% Seq Cov

N	Q	S	A	L	T	Q	P	R	S	V	S	G	S	P	G	Q	S	V	T	I	S	C	T	G	T	25	
26	S	S	D	I	G	G	Y	N	F	V	S	W	Y	Q	Q	H	P	G	K	A	P	K	L	M	I	50	
51	Y	D	A	T	K	R	P	S	G	V	P	D	R	F	S	G	S	K	S	G	N	T	A	S	L	75	
76	T	I	S	G	L	Q	A	E	D	E	A	D	Y	Y	C	C	S	Y	A	G	D	Y	T	P	G	100	
101	V	V	F	G	G	G	T	K	L	T	V	L	G	Q	P	K	A	A	P	S	V	T	L	F	P	125	
126	P	S	S	E	E	L	Q	A	N	K	A	T	L	V	C	L	L	I	S	D	F	Y	P	G	A	V	150
151	T	V	A	W	K	A	D	S	S	P	V	K	A	G	V	E	T	T	T	P	S	K	Q	S	N	175	
176	N	K	Y	A	A	S	S	Y	L	S	L	L	T	P	E	Q	W	K	S	H	R	S	Y	S	C	Q	200
201	V	T	H	E	G	S	T	V	E	K	T	V	A	P	T	E	C	S	C								


Light Chain Combined Coverage

70% Sequence Coverage


- ETD 15ms, 25ms
- HCD 10
- CID 30

N **Q** S A L T Q P R S V S G S P G Q S V T I S C T G T 25
26 S S D I G G Y N F V S W Y Q Q H P G K A P K L M I 50
51 Y D A T K R P S G V P D R F S G S K S G N T A S L 75
76 T I S G L Q A E D E A D Y Y C C S Y A G D Y T P G 100
101 V V F G G G T K L T V L G Q P K A A P S V T L F P 125
126 P S S E E L L Q A N K A T L V C L L I S D F Y P G A V 150
151 T V A W K A D S S P V K A G V E T T T P S K Q S N 175
176 N K Y A A S S Y L L S L L T P E Q W K S H R S Y S C Q 200
201 V T H E G S T V E K T V A P T E C S C

Heavy Chain: ETD HD Fragmentation (CS 52+)



Heavy Chain HCD Fragmentation (CS 52+)



Combined Sequence Coverage

31% Sequence Coverage

```
N  E V Q L V E S G G G L V Q P G G S L R L S C V A S 25
26  G F T L N N Y D M H W V R Q G I G K G L E W V S K 50
51  I G T A G D R Y Y A G S V K G R F T I S R E N A K 75
76  D S L Y L Q M N S L R V G D A A V Y Y C A R G A G 100
101 R W A P L G A F D I W G Q G T M V T V S S A S T K 125
126 G P S V F P L A P S S K S T S G G T A A L G C L V 150
151 K D Y F P E P V T V S W N S G A L T S G V H T F P 175
176 A V L Q S S G L Y S L S S V V T V P S S S L G T Q 200
201 T Y I C N V N H K P S N T K V D K K V E P K S C D 225
226 K T H T C P P C P A P E L L G G P S V F L F P P K 250
251 P K D T L M I S R T P E V T C V V V D V S H E D P 275
276 E V K F N W Y V D G V E V H N A K T K P R E E Q Y 300
301 N S T Y R V V S V L T V L H Q D W L N G K E Y K C 325
326 K V S N K A L P A P I E K T I S K A K G Q P R E P 350
351 Q V Y T L P P S R D E L T K N Q V S L T C L V K G 375
376 F Y P S D I A V E W E S N G Q P E N N Y K T T P P 400
401 V L D S D G S F F L Y S K L T V D K S R W Q Q G N 425
426 V F S C S V M H E A L H N H Y T Q K S L S L S P G C
```

- ETD 15ms, 25ms
- HCD 10%