

How to provide a saliva sample

Capture the QR code on the bag and log into <https://rtr.ipo.rutgers.edu/questionnaire> to claim your test kit

STEP 1

- Complete the online questionnaire and print the requisition. Insurance information must be entered.
- Labels (3 types): where to apply
 - ✓ On bag: netid (write in) AND white test kit barcode
 - ✓ On requisition: test kit barcode
 - ✓ On tube: name/date of birth/netid (see step 5)



STEP 2


Fill the tube with saliva to the wavy black line (not including bubbles). Do NOT overfill.

STEP 3

Replace the funnel with the fluid cap and screw tightly until the blue fluid releases into the tube.

STEP 4

Shake the tube for at least 5 seconds. Make sure the blue fluid releases.



STEP 5

Write your name, date of birth, and netid on the blank label and apply to the tube. Do NOT cover the test tube barcode.



STEP 6

Place the tube in the bag and the requisition form in the outside pocket.

Do NOT include the funnel or packaging.



Please produce and submit the sample in one day. Do not eat, drink or smoke for 30 minutes before producing the sample.

Please drop off your sample in a COVID test drop off box, Sunday-Thursday, last weekly pickup Friday 8 am.